


COALITION™ ACDF System

The *MILDER*™* ACDF



DTS Guide for One Step Insertion



Adjacent Segment Treatment



* Minimally Invasive, Less Disruptive, Earlier Recovery

COALITION™ ACDF System Benefits

COALITION™ Spacer

- Anatomical profile plate, screw and spacer system
- Three axial footprints (12 x 14mm, 14 x 16mm, 15 x 18mm)
- Two sagittal profiles (parallel and 7° lordotic)
- Seven heights 6mm-12mm
- Torsional stabilization features
- Pre-assembled locking mechanism in plate
- Audible, tactile and visual confirmation of screw blocking



Variety of Screw Sizes


- 3.6mm or 4.2mm diameter
- Self-drilling or self-tapping
- Lengths from 12mm to 20mm
- Variable and fixed angle screws



Screw Angulation

- Screw orientation is 35° cephalad/caudal
- Variable angle screws up to $\pm 4^\circ$
- Medial orientation: 12°-16°

Variable Angle Screw Fixed Angle Screw



COALITION™ Instrumentation

- Minimally invasive system and technique
- Streamlined holder for minimal access
- Pre-set angle DTS guide holder allows one step insertion
- Contained CV joint and straight instruments for ease of insertion



www.globusmedical.com

Globus Medical
Valley Forge Business Center
2560 General Armistead Avenue
Audubon, PA 19403

CE REP: RMS - UK Limited
28 Trinity Road, Nailsea, Somerset, BS484NU England

CE 0297

Customer Service:
Phone 1-866-GLOBUS1 (1-866-456-2871)
Fax 1-866-GLOBUS3 (1-866-456-2873)

©2009 Globus Medical. All rights reserved. COALITION is a trademark of Globus Medical. Please refer to the product insert for description, indications, contraindications, warnings, precautions, and other important information.

GMSPS55
08.09