

REVOLVE[®] Pedicle Screw System

Multi-Level Capability



Integrated Rod Reduction



Dependable
Screw
Connection



Non-Threaded Locking Cap




REVOLVE[®] System Features

The easiest pedicle screw system—now in MIS.

Every step of the minimally invasive procedure has been enhanced.

REVOLVE[®] System

- Multi-level capability
- Powerful rod reduction up to 22mm
- Strong connection between screw and reduction sleeve
- Non-threaded locking cap
- Easy to use instrumentation
- Low profile implants and instruments



Polyaxial Screws/Cannulated Screws

- $\pm 30^\circ$ screw angulation (60° total) provides intra-operative flexibility
- Low profile top-loading screw design
- Allows for compression and distraction while screw angle is locked
- Low torque locking mechanism
- Self-tapping design
- Screw diameters 5.5mm, 6.5mm, 7.5mm, 8.5mm
- Screw lengths from 30mm to 90mm



Rods

- 5.5mm rod diameter
- Rod lengths from 80mm to 300mm
- Straight rods can be intra-operatively contoured



Reduction Sleeve

- Integrated threads allow for up to 22mm in reduction
- Large slots along the sides enable easy targeting for rod entry
- Threaded reducer powerfully reduces rod
- Cap insertion with only 90° finger turn



www.globusmedical.com

Globus Medical
Valley Forge Business Center
2560 General Armistead Avenue
Audubon, PA 19403

CE REP: RMS – UK Limited
28 Trinity Road, Nailsea, Somerset, BS484NU England

CE 0297

Customer Service:

Phone 1-866-GLOBUS1 (1-866-456-2871)

Fax 1-866-GLOBUS3 (1-866-456-2873)

©2009 Globus Medical. All rights reserved. REVOLVE is a registered trademark of Globus Medical. Please refer to product insert for description, indications, contraindications, warnings, precautions, and other important information.

GMSP53
11.09