

SIGNATURE™ TLIF System

A PATRIOT® Spacer

Unprecedented Control



MIS TLIF Compatible




Easy Insertion



SIGNATURE™ TLIF System *A PATRIOT® Spacer*

SIGNATURE™ Spacer


- Integrated pivoting mechanism for self-articulation
- Two axial footprints 10x28mm and 11x33mm
- Heights of 7–13mm, 15mm and 17mm
- Easy insertion
- Curved geometry to match vertebral anatomy
- Large axial and anterior holes for bone ingrowth
- Tantalum radiographic markers help navigate ideal cage placement



SIGNATURE™ Instruments

- Inserter:
 - Low profile design for minimally invasive surgery
 - Maintains implant connection through final placement for optimal control
 - Lateral stabilizers engage implant for insertion and release for rotation
- Paddle distractors, scrapers and tamps facilitate discectomy and insertion
- Trials and rasps mimic implant geometry

Height	Small	Large
7mm	368.207	368.307
8mm	368.208	368.308
9mm	368.209	368.309
10mm	368.210	368.310
11mm	368.211	368.311
12mm	368.212	368.312
13mm	368.213	368.313
15mm	368.215	368.315
17mm	368.217	368.317



www.globusmedical.com

Globus Medical
Valley Forge Business Center
2560 General Armistead Avenue
Audubon, PA 19403

Customer Service:

Phone 1-866-GLOBUS1 (1-866-456-2871)

Fax 1-866-GLOBUS3 (1-866-456-2873)

CE REP: RMS - UK Limited
28 Trinity Road, Nailsea, Somerset, BS484NU England



©2009 Globus Medical. All rights reserved. SIGNATURE is a trademark of Globus Medical.

Please refer to the product insert for description, indications, warnings, precautions, and other important information.

GMSPS37
03.09