

SUSTAIN[®]-0


Teeth to Resist Expulsion




Large Axial Hole



Self-Distraction



Tapered Leading Edge





SUSTAIN®-0 Features

- Tapered leading edge for easier insertion
- Rounded corners allow for rotation during insertion
- Grip-style connection for controlled implant positioning
- Teeth to resist expulsion
- 1mm height increments for accurate fit
- 22mm, 26mm and 30mm lengths for better endplate coverage
- 8mm, 10mm and 12mm width options for an ideal fit

SUSTAIN®-0 Small Implant Set 904.903* SUSTAIN®-0 Implant Set 904.908

Height	8x22mm*	10x22mm*	10x26mm	10x30mm	12x26mm	12x30mm
8mm	304.078	304.088	304.008	304.808	304.968	304.868
9mm	304.079	304.089	304.009	304.809	304.969	304.869
10mm	304.080	304.090	304.010	304.810	304.970	304.870
11mm	304.081	304.091	304.011	304.811	304.971	304.871
12mm	304.082	304.092	304.012	304.812	304.972	304.872
13mm	304.083	304.093	304.013	304.813	304.973	304.873
15mm	304.085	304.095	304.015	304.815	304.975	304.875
17mm	304.087	304.097	304.017	304.817	304.977	304.877

*Implants highlighted in gray are additionally available.



Indications: The SUSTAIN® Spacer is a vertebral body replacement device intended for use in the thoracolumbar spine (T1-L5) to replace a collapsed, damaged, or unstable vertebral body due to tumor or trauma (i.e., fracture). The SUSTAIN® Spacer is intended to be used with supplemental spinal fixation systems that have been labeled for use in the thoracic and/or lumbar spine (i.e., posterior pedicle screw and rod systems, anterior plate systems, and anterior screw and rod systems). The interior of the spacer can be packed with bone grafting material. The SUSTAIN® Spacer is designed to provide anterior spinal column support even in the absence of fusion for a prolonged period.



www.globusmedical.com

Globus Medical
Valley Forge Business Center
2560 General Armistead Avenue
Audubon, PA 19403

Customer Service:
Phone 1-866-GLOBUS1 (1-866-456-2871)
Fax 1-866-GLOBUS3 (1-866-456-2873)

CE 0297

©2009 Globus Medical. All rights reserved. SUSTAIN is a registered trademark of Globus Medical. Please refer to product insert for description, indications, contraindications, warnings, precautions, and other important information.

GMSPS18
03.09